

CLEANING & DISINFECTION OF SUBMERSIBLE INSTRUMENTS

Pro-Cleanse™

Triple Action Enzymatic Cleanser




- Ideal for pre-cleaning of soiled equipment in preparation for disinfection and sterilisation
- Effective on dried deposits of blood and bodily fluids
- Suitable for use on stainless steel, metal, glass and plastics
- Environmentally friendly and readily biodegradable

Cleaning



1. Fill a tray (or other suitable container) with Vet-Cleanse Enzymatic Cleanser diluted at a concentration of 3-7.5ml per 1Ltr of water



2. Submerge soiled instruments entirely in the tray and soak for 20 minutes



3. Remove the instruments from the tray and rinse them thoroughly with purified water



4. Clean the instruments with a high level disinfectant solution such as Pro-Cleanse Instrument Disinfectant or using an autoclave

Pro-Cleanse™

Instrument Disinfectant



- Broad spectrum biocide highly effective against bacteria, viruses, fungi and yeasts
- Cleans and protects without abrasion, bleaching, odour or taint
- Suitable for use on all materials including metal, glass, plastics, rubbers and polymers

Disinfecting



1. Pre-clean instruments using Pre-Cleanse Enzymatic Cleanser and rinse thoroughly



2. Fill a tray (or other suitable container) with Pro-Cleanse Instrument Disinfectant diluted at a concentration of 1:100 with water



3. Submerge the instruments entirely in the tray and soak for 10 minutes



4. Remove the instruments from the tray and rinse them thoroughly with purified water